



PTT-B54

Injection in a Box™

Patent Nos. US10,217,380, US10,535,282, US10,586,470, US11,195,434,
US11,195,435, Patents Pending

TECHLINE® Device Specifications

March 11, 2022

The Injection in a Box™ task trainer is for training students in intramuscular injection. The simulated SOFTTECH® skin provides realistic feedback to the trainee, while the injected fluid is collected by the removable sponge.

Construction-

The representative skin and tissue of the PTT-B54 is constructed of high-quality silicone rubber, with layers of coloration that provide a depth and realism to the product. In addition, the SOFTTECH® material has a realistic feel, and is very durable and UV resistant.

Operation-

The Injection in a Box™ can be used for administration of intermuscular injections. Can inject fluid to simulate IV administration. The moisture collecting sponge allows for the collection of the injected fluid in the cavity. Can be used inside or outside of the case.

Maintenance-

The device should be rinsed in clean, warm water and allowed to dry before storage. The simulation may be removed from the storage case for cleaning. The device may be washed with mild detergents like dish soap or cleaned with a disinfecting agent like isopropyl alcohol or a diluted bleach solution, as it is durable, and chemically resistant. Allow the device to dry fully before storage to prevent mold or mildew growth.



PTT-B54 Injection In A Box™